



TOTAL ORTHOPAEDIC CARE

CONTACT US:

****** +254 703 776 694

www.orthomedics.co.ke



1st Floor Room 102 Sagaas Centre, Along Makasembo Road, Opp. Kmtc Adjacent To Aic Fellowship. P.O Box 3705 - 30100 Eldoret.

TABLE OF CONTENTS

ACHILLIES BOOT	
ARTHROSCOPY KNEE SUPPORT	
PATELLA AND LIGAMENT KNEE SUPPORT LONG OPEN	
PATELLA AND LIGAMENT SUPPORTED KNEE SUPPORT	
R.O.M WALKER	
CLOSED BASIC KNEE SUPPORT	
STEEL HINGE KNEE SUPPORT OPEN FRONT	
MALLEOL SUPPORTED ANKLE SUPPORT	
HILUX VALGUS NIGHT SPLINT SOFT	
REFLEKS AFO RIGHT & LEFT	
UNIV MALLEOLAR ANKLE SUPPORT	
THIGH LENGTH ANTIEMBOLIC STOCKING	
THUMB SUPPORTED WRIST SPLINT RIGHT & LEFT	
SILICONE BREAST PROSTHESIS	
LUMBOSACRAL CORSET 26CM	



OL - 1104 Achilles Boot

Beden/Size	Ölçüler (Num.)
	34 -37
	38 - 41
L	42 - 45



Internal surface is made of towel foam fabric. Is the ankle orthesis helping early mobilization and ankle plantar flexion height can be adjusted by 3-stage pads in heel



OL - 2108 Arthroscopy Knee Support

Beden/Size	Ölçüler (cm)
S	31 x 34
М	34 x 37
L	37 x 40
XL	40 x 44
XXL	44 x 48



Is made of neoprene material and designed anatomically. Pocket system inside is designed for gel used in hot and cold application. Patient being operated and old can use it easily since it is opened fully.

INDICATIONS: Used in the case that hot and cold application is preferred for treatment of post-arthroscopy edema and pain



OL - 2103 A
Patella and
Ligament
Supported knee
support
(Long) (Open)

Beden/Size	Ölçüler (cm)
S	31 x 34
М	34 x 37
L	37 x 40
XL	40 x 44
XXL	44 x 48



Is made of neoprene material, designed anatomically. Adjacent flexible underwire increases stabilization and relieve the burden on the knee.

Allows for easy use in upper calf by means of front velcro. Circumference of patella is open and supported with a pad. Use of it is easy since the part of upper calf can be opened.

INDICATIONS: Used in the cases of arthritis, intra-articular fluid increases, patella pains, for post-operation protection and rehabilitation, heat and compression.



OL - 2102 Patella and ligament supported knee support

Beden/Size	Ölçüler (cm)
S	31 x 34
М	34 x 37
L	37 x 40
XL	40 x 44
XXL	44 x 48



Is made of neoprene material, designed anatomically. Adjacent flexible underwire increases stabilization and relieve the burden on the knee.

Allows for easy use in upper calf by means of front velcro. Circumference of patella is open and supported with a pad.

INDICATIONS: Used in the cases of arthritis, bursitis, intra-articular fluid increases, patella pains, meniscus, first degree sprain and ligament instabilities, for post-operation protection.



OL - 1105 R.O.M. Walker

Beden/Size	Ölçüler (Num.)
	34 -37
	38 - 41
	42 - 45



Inside of it is made of cotton fabric and helps early mobilization as ankle dorsi flexion and plantar flexion degree can be adjusted.

INDICATIONS: Tendon-ligament soft tissue injuries, chronic ligament instability, distal tibia and fibula fractures not being instable or with internal fixator, Achilles Tendon rupture, stabile foot and ankle fractures, for post-op stabilization.



OL-2100 Closed basic knee support

	Ölçüler (cm)
S	31 x 34
М	34 x 37
	37 x 40
	40 x 44
XXL	44 x 48



Is made of neoprene material and designed anatomically. Allows for easy use in upper calf by means of front velcro. It supports knee joint. It reduces the pain by relieving the burden on knee.

INDICATIONS: used in the cases of arthritis, intra articular fluid increase, for post-operation protection and rehabilitation, heat and compression.



OL-2104 **Steel Hinged** knee support

Beden/Size	Ölçüler (cm)
S	31 x 34
М	34 x 37
L	37 x 40
XL	40 x 44
XXL	44 x 48



Is made of neoprene material and designed anatomically. Steels on both sides provide maximum support and allow for controlled motion. There is velcro increasing the stabilization in upper and lower calf.

INDICATIONS: Used in the cases of arthritis, bursitis, tendinitis, sprain and strain, meniscus problems, ligament instability, post-sport injury and post-op protection.



OL - 2403 **Malleol Supported** Ankle Support

,	Jupp
Beden/Size	Ölçüler (cm)
S	21 x 23
М	23 x 25
L	25 x 27
XL	27 x 29
XXL	29 x 31



Is made of neoprene material. Internal and external malleolar are supported with anatomically designed pads. It supports ankle and reduces edema and

INDICATIONS: Mild ligament injuries and instability, joint inflammation and arthritis, for post-op and injury protection and support.



OL - 30 Hallux Valgus Night Splint (Soft)

STANDARD SIZE

RIGHT AND LEFT IS

AVAILABLE

Inside of it is made of towel fabric and air-permeable material. It prevents deviation of toe to lateral. It is used for Hallux Valgus deformity especially at night.



Is made of thermoplastic material and internal surface is coated with plastozote. It keeps ankle at 90° and velcro around ankle offers easy use. Can be used in the shoes.

INDICATIONS: Used in the cases that ankle should be at 90°, low foot deformity in the event of nerve injuries such as hemiplegia, herniated disk



STANDARD SIZE

Unv. Malleolar ankle support



OL - 107 STANDARD.

Seating Cushion



Knee length compression stockings



Compression stockings thigh length



Thigh length antiemboli stockings



Is made of neoprene material. Its use is easy by means of the zip.

INDICATIONS: Used for distension and inguen problems, slimming of waist, hip and thigh and for heat and compression.



OL - 7011 STANDARD SIZE

Unv. Thumb supported wrist splint



OL - 21 Unv. Wrist splint with Thumb

STANDARD SIZE.

RIGHT AND LEFT IS

AVAILABLE

Is made of velcro neoprene fabric and circumference of wrist is reinforced with flexible velcro.

INDICATIONS: Used in the cases of Dequarvain tenosynovitis, carpal tunnel and Dequarvain, primary CMC and MCP joint degenerations, for immobilization of aching wrist and thumb, for protection after injury and operation, and in the post-plaster process for fractures.



OL - 7010 Unv. Wrist splint

STANDARD SIZE



OL - 15 N Neoprene Dequervain Splint

Beden/Size	Ölçüler (cm)
S	14 x 16
М	16 x 18
L	18 x 20
XL	20 x 22
XXL	22 x 24

RIGHT AND LEFT ARM IS AVAILABLE.



Is made of neoprene material and has an aluminum support restricting the thumb motion with adjustable slope.

INDICATIONS: Used in the cases of Dequarvain syndrome, primary CMC and MCP joint irritations, for protection for thumb ligament laxity, heating and compression.



OL - 300 Stimson Bandage

n 1 (n)	21 m / 3
	Ölçüler (cm)
S	34 x 38
	38 x 44
	44 x 50
	50 x 58
XXL	58 x 66



Internal surface of which is made of towel and foam fabric. Should pad ensures uniform distribution of the burden.

INDICATIONS: used in the cases of Achromia, clavicular, dislocation-semi dislocation.



OL - 21 Thumb supported Wrist splint

Beden/Size	Ölçüler (cm)
S	14 x 16
М	16 x 18
L	18 x 20
XL	20 x 22
XXL	22 x 24

RIGHT AND LEFT ARM IS AVAILABLE.



Is made of antiperspirant cotton fabric. It has an aluminum support restricting motions of wrist and thumb. There is a bandage providing stabilization around the wrist.

INDICATIONS: Used in the cases of Dequarvain tenosynovitis, carpal tunnel and Dequarvain, primary CMC and MCP joint degenerations, for resting in cases that wrist and thumb is aching, for protection after operation and injury



OL - 19 Long Wrist Splint

Beden/Size	Ölçüler (cm)
S	14 x 16
	16 x 18
	18 x 20
	20 x 22
XXL	22 x 24

RIGHT AND LEFT ARM IS AVAILABLE.



Is made of antiperspirant cotton fabric. Its aluminum support supporting wrist from volar face restricts wrist motions. Fingers are free from MCP joint. There is a bandage providing stabilization around the wrist.

INDICATIONS: Used in the cases of carpal tunnel syndrome, tendinitis, tenosynovitis, tenovaginitis, for resting in the case of wrist pain, protection after injury and operation.



OL - 19 Long Wrist Splint

Beden/Size	Ölçüler (cm)
S	14 x 16
М	16 x 18
L	18 x 20
XL	20 x 22
VVI	22 v 24

RIGHT AND LEFT ARM IS AVAILABLE.



Is made of antiperspirant cotton fabric. Its aluminum support supporting wrist from volar face restricts wrist motions. Fingers are free from MCP joint. There is a bandage providing stabilization around the wrist.

INDICATIONS: Used in the cases of carpal tunnel syndrome, tendinitis, tenosynovitis, tenovaginitis, for resting in the case of wrist pain, protection after injury and operation.



OL - 2009 Lymph edema arm support

Used for edema cases of arm. A model starting from the wrist to the thumb is available. By means of high air-permeable feature, it maintains body temperature. Elastic and soft weaving allow for easy use.



Oval silicone exercise balls sued for exercise in rehabilitation. 3 hardness levels are available as red (hard), blue (medium), green (soft).



OL - 35 Dynamic splint

Beden/Size	Ölçüler (cm)
S	16 x 18
	18 x 21
	21 x 24

RIGHT AND LEFT IS AVAILABLE



OL - 96 Silicone Exercise ball


Is made of thermoplastic material and is a hand and wrist orthesis internal surface of which is coated with towel fabric. It allows for active finger flexion. It keeps wrist and fingers at extension passively.



OL - 3020 Silicone Breast prosthesis

USED BY 15 HEIGHTS (60-130) SIZES.

Is made of high-quality silicone and used when breast tissue is excised. It gives a natural appearance and reduces pains due to muscle imbalance.



OL - 3020 Breast prosthesis bra

USED BY15 HEIGHTS (60-130) SIZES.

Is produced for breast prosthesis. It supports the prosthesis and weight and allows comfortable use and gives confidence. Internal surface is made of non-irritant cotton fabric. Shoulder slings are softened for easy use. The prosthesis can be used in both directions by means of pockets on both sides.



OL - 320 K Lumbosacral Corset (32 cm) With belt

Beden/Size	Ölçüler (cm)
S	75 x 85
М	85 x 95
L	95 x 105
XL	105 x 115
XXL	115 x 125

Is made of antiperspirant air-permeable fabric and has steel and flexible underwire. Cross band extra belt support relieves the burden on disks. It is 32 cm and supports the whole waist and lower back.

INDICATIONS: Used in the cases of lumbago, acute or chronic lumbar pains, degenerative disk lesions, spondylosis, for post-op immobilization.



OL - 260 K Lumbostad Corset (26 cm) With belt

Beden/Size	Ölçüler (cm)
S	75 x 85
М	85 x 95
	95 x 105
	105 x 115
XXL	115 x 125
1	



Is made of antiperspirant air-permeable fabric and has steel and flexible underwire. Cross band extra belt support relieves the burden on disks.

INDICATIONS: Used in the cases of lumbago, acute or chronic lumbar pains, degenerative disk lesions, spondylosis, spondylolisthesis, for post-op mobilization.



OL - 2801 STANDARD SIZE

Shoulder support



OL - 1001 Hyper extension corset

S	40 x 49
М	50 x 60



Is made of special aluminum material and surfaces contacting the skin is coated with towel pad. Based on the principle of 3-point, it restricts flexion of the body and stimulates hyper extension position. It makes spine stabilize against lateral flexion and rotation motions. By means of the screw system on it, length and width can be adjusted.

INDICATIONS: Used in the cases of spine fractures and traumas, osteoporosis, compression fractures, kyphosis and post-op process.



OL - 3005 Forearm crutch



OL - 3004 P **Underarm crutch**

Is made of aluminum pipe, length of stick is adjustable and has a vacuum rubber tire. Hand holding and arm supporting part designed anatomically is one-component and is made of high-quality rigid plastic.

Is made of aluminum pipe, length of stick is adjustable. The part under arm and hand holding part are made of PVC. Since there is a vacuum rubber tire under the crutch, its stabilization is quite high. 3 lengths are available



OL - 3000 Luks Mobile Walker


Is made of aluminum pipe and heights of four legs are adjustable. By means of joints, it can be moved back and forth while walking. It can be folded when not use and stored with no space



SABO Classic For Man

\$SMG1391 OSMG1392





Is made of silicone material coated with fabric. It reduces burden on the heel causing pain.

INDICATIONS: Pes planus (flat foot) Pes cavus (supinated foot Used for supporting longitudinal arch.



OL - 100 Lüks Silicone Lux Insole

Beden/Size	Ölçüler (cm)
XS	33 x 34
	35 x 36
М	37 x 38
	39 x 40
XL	41 x 42
XXI	43 x 44



Is made of silicon material and provides balanced distribution of body burden. Absorbs ground reaction force and reduces burden and stress on joints. Arch is elevated completely.

INDICATIONS: Flat-footed, plantar fasciitis, heel pains, pain of metatarsus



Meets the requirement of shoes when the leg is encased in plaster.



OL - 3009 Lux Folding toilet seat riser



OL - 3009 Folding toilet seat riser

Meets the need to use the toilet of patients having walking difficulty at closed areas such as home, hospital. It is practical since the patient does not need to be lifted and it can be removed without seeing potty. It can be stored by folding when not used

Meets the need to use the toilet of patients having walking difficulty at closed areas such as home, hospital. It is practical since the patient does not need to be lifted and it can be removed without seeing potty. It can be stored by folding when not used



OL - 3010 Wheelchair



OL - 3010 E Domestic Wheelchair

It can be folded and foot switches can be removed, is chrome plated, side armrests can be removed without lifting the patient and is available in filled and air wheel models.

Is the wheelchair with potty that can be used when the patient cannot stand in closed areas such as house, hospital.



OL - 10

Portable Toilet Seat riser (Cover)



REGD. OFFICE:

1ST FLOOR ROOM 102 SAGAAS CENTRE,
ALONG MAKASEMBO ROAD,
OPP. KMTC ADJACENT TO AIC FELLOWSHIP.
P.O BOX 3705 - 30100 ELDORET.

FOR INQUIRY: +254 703 776 694